

HairMax LaserComb Ultima 9 Classic & Ultima 12

Quick guide

The HairMax LaserComb® is a very effective device when used correctly. It is manufactured using precision components and should last 10-15 years under normal usage. Please take time to read the entire User Manual before operating the LaserComb.

HairMax LaserComb Ultima 9 Classic & Ultima 12

Use the LaserComb Ultima 9 Classic for 11 minutes, a LaserComb Ultima 12 for 8 minutes 3 times a week. The treatments should be spaced every other day, for example Monday, Wednesday, Friday. Move the LaserComb slowly 1/2 inch every 4 seconds.

It is recommended, that the LaserComb be used on a clean scalp after washing. Your hair can be wet, towel dried or dry. The LaserComb should be used prior to applying any topical treatments, gels, sprays or concealing products (they can hinder the flow of laser light).

- Charge your LaserComb by connecting the charger to the stand (figure A) or by placing the power cord directly into the handle end of the LaserComb (Figure B). For initial charging, allow 2-3 hours to fully charge your LaserComb. The indicator light on the power button will flash red during the charging process. Once charging is complete, the light will stop flashing and remain on.
- Turn on the LaserComb by pressing the power button and hold until power button light turns on. Next, press the laser button until it lights up. That means the LaserComb is ready for use(Figure C).

Laser Light - Avoid direct eye exposure.

- Before beginning treatment, be sure to test the LaserComb for sensitivity (see page 5 of your User Manual). To start your treatment place the LaserComb flat on your scalp so that both rows of teeth are touching the scalp (Figure D).
- Start at your hairline. Move the LaserComb from front to back of scalp in a slow combing motion. Stop every ½ inch (1cm) and pausing for 4 seconds, until you hear a vibration signal (Figure E). Do not lift the LaserComb off your scalp while combing from front to back.
- Lift the LaserComb off the scalp, then start the treatment on the side of the head, pointing the comb from the bottom up (Figure F). Continue this operation from front to back, then side to side, covering the entire head, respecting the 4-second rule and gentle movements. When the LaserComb is on it will beep and vibrate every four (4) seconds. This means that it is time to move the LaserComb to the next area of the scalp. Repeat these passes for the total treatment of 8 minutes for Ultima 12 and 11 minutes for Ultima 9 Classic.
- If you have long hair, use your fingers to lift your hair out of the way while using the LaserComb to be sure that the laser energy reaches your scalp (Figure G).
- After treatment, turn off the LaserComb by pressing and holding the power button (lower button) until the indicator light turns off(Figure H). Store the device in a cover or leave it in

a stand. You can connect the power cable which will allow the laser to work on the next treatment.

The HairMax LaserComb is indicated to treat Androgenetic Alopecia, and promote hair growth in males who have Norwood Hamilton Classifications of IIa to V and in females who have Ludwig (Savin) I-4, II-1, II-2, or frontal patterns of hair loss and who both have Fitzpatrick Skin Types I to IV.

Please Note:

- The time required for treatment with HairMax LaserComb is 8 minutes for Ultima 12 and 11 minutes for Ultima 9 Classic.
- It is important to move the LaserComb very gently. If the laser is moved too quickly, the treatment may not produce optimal results.
- The LaserComb will turn off automatically after 5 minutes.
- Although there are no negative effects of abuse, it is recommended to use the LaserComb according to the instructions.
- LaserComb has advanced electronic circuits to ensure that the device works properly. If for any reason you have problems with the device, before sending it to us for repair, please contact us by phone or email.