

# HairMax LaserBand 41

## Quick User Guide

Use your HairMax for as little as 3 minutes, 3 times a week, ideally every other day (for example, Monday, Wednesday, Friday).

### FOR MAXIMUM BENEFIT:

- scalp should be freshly washed (hair can be wet, dry or damp)
- use your HairMax prior to applying styling products (these can build up and impede the laser light)
- regularly clean the laser area with a damp cloth. do not use harsh detergents or cleaning products on any part of your HairMax

### CHARGING OPTIONS

#### Charge your HairMax for 4-5 hours before initial use.

Prior to charging your HairMax, make sure your HairMax is turned off.

**OPTION 1** Plug the power cord into the charging stand and the other end into an electrical wall outlet (Figure A). Place your HairMax on the charging stand. And wait until the indicator light starts flashing. A steady light indicates your HairMax is fully charged.

**OPTION 2** Place the power cord directly into the charging socket on the bottom left side of your HairMax and the other end into an electrical wall outlet (Figure B).

And wait until the indicator light starts flashing. A steady light indicates your HairMax is fully charged (Figure C).

**Do not use your HairMax while it is plugged into the power cord.**

Prior to treatment, test your HairMax for sensitivity (Figure D). For more information, see page 5 of the User Manual.

### TREATMENT METHOD

Your HairMax will be used in 6 positions on your scalp. Leave it in place for a minimum of 30 seconds for each position. It will vibrate every 30 seconds to let you know when to move it to the next position.

To turn on your HairMax, press and hold the power button until the red lasers turn on (Figure E). Next, place it on your frontal hair line. This is position 1. When you feel the vibration, slide it backwards about 2 inches (5 cm) to position 2. Leave it in position 2 until you feel the vibration again. Repeat this process for all six positions on your scalp. The entire treatment should be a minimum of 3 minutes.

After treatment, Turn off your HairMax by pressing and holding the power button until the lasers turn off (Figure E). Remove your HairMax from your head and place it on the charging stand for future use.

**Please Note:** Your HairMax turns off automatically after 12 minutes.

The HairMax® LaserBand is a highly effective device when used correctly. It is manufactured using precision components and should last 10-15 years under normal usage. Please take time to read the User Manual before operating your HairMax.

Laser Light - Avoid direct eye exposure.

The HairMax LaserBand is indicated to treat Androgenetic Alopecia, and promote hair growth in males who have Norwood Hamilton Classifications of IIa to V and in females who have Ludwig (Savin) I-4, II-1, II-2, or frontal patterns of hair loss and who both have Fitzpatrick Skin Types I to IV.